

Hygicult Y&amp;F (Easicult M)

Hygicult E/ $\beta$ -GUR

Hygicult E

Hygicult TPC

# Simple surface hygiene monitoring.

## Hygicult® Model Chart

## Hygicult®

For monitoring microbial hygiene in solid and liquid materials.

Storage at 18–25°C avoiding draughts and temperature fluctuations.

Detection range for bacteria is  $10^3$ – $10^7$  CFU/ml in liquids or 1–100 CFU/cm<sup>2</sup> on surfaces.

### Hygicult TPC (Cat. no. 68010)

- Shelf life of six months from the date of manufacture.
- Incubation at 35–37°C for 24 h, at 27–30°C for 48 h, or 18–25°C for up to five days.

### Hygicult E (Cat. no. 68012), Hygicult E/ $\beta$ -GUR (Cat. no. 68267), Hygicult CF (Cat. no. 06036)

- Shelf life of five months from the date of manufacture.
- Incubation at 35–37°C for 24–48 hours.

### Hygicult Y&F (Easicult M) (Cat. no. 68013)

- Shelf life of six months from the date of manufacture.
- Incubation at 27–30°C for three to five days.
- Detection range for yeasts is  $10^3$ – $10^6$  CFU/ml in liquids or 1–80 CFU/cm<sup>2</sup> on surfaces.

For more detailed information, please consult the instructions for use.

It's so easy to use!



Dip, press or swab



Incubate



Interpret


# Hygicult®

## Colony densities




### Hygicult TPC

Total bacteria,  
yeasts & moulds




### Hygicult E Hygicult CF

Total bacteria,  
yeasts & moulds



### Hygicult E/β-GUR

Growth on β-GUR medium



### Hygicult Y&F (Easicult M)

Yeasts & Moulds

